# C-Bond®

- Seite 2 - de: Gebrauchsanweisung  
Page 3 - en: Instructions for use  
Page 4 - fr: Manuel d'instruction  
Página 5 - es: Instrucciones de uso  
Pagina 6 - it: Istruzioni per l'uso  
Strana 7 - bs: Uputstvo za upotrebu  
Pagină 8 - ro: Instrucțiuni de utilizare  
Sayfa 9 - tr: Kullanma kılavuzu  
Strona 10 - pl: Instrukcja użytkowania  
Стр. 11 - ru: Инструкция по применению

*We know how*

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**Patientenzielgruppe**

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## Instructions for use

**Patient target group**

Children, adolescents, adults and older patients

**Intended users**

Dentist

**Composition**

Qualitative data:

Oligomeric methacrylates, HEMA, camphorquinone

Quantitative data in % G/G: Methacrylates 98, Additive 2

**Intended purpose of the product**

C-Bond is an adhesive with excellent adhesion to dentin and etched enamel. It improves the adhesion between composite and tooth surface.

**Indications**

C-Bond is used for improved adhesion of composites to enamel and dentin. It also improves the adhesion between glass ionomer underfilling cements and composite fillings ("sandwich technique").

**Contraindications**

C-Bond must not be used in the pulp area.

**Performance characteristics of the product**

Curing by irradiation with light.

**Specific precautions**

None

**Application**

The cavity should be prepared according to the rules of the adhesive filling technique. Depending on the situation, calcium hydroxide paste (Calcident 450) or an underfilling cement should be placed in deep cavities. According to the etch-bond technique, the enamel is etched, rinsed and then dried. Glass ionomer cement filling should not be etched, but only roughened.

After the prepared cavity is dry, 1-2 drops of light-curing C-Bond are dropped onto the mixing block (reseal the bottle immediately) and, using a brush or foam pellets, laid down in a thin layer. Excess is removed with a dry brush or a pellet. This is followed by curing with halogen light (15 seconds).

Light devices should emit at 450 nm at 500 mW/cm<sup>2</sup> and should be checked regularly.

**Warnings****Residual risks**

The processing instructions and precautionary measures described in the instructions for use of the product must be observed in any case. If, despite following these instructions, an adverse reaction occurs in the patient or other adverse events occur, we request that you inform us in every case, describing the accompanying circumstances as precisely as possible.

**Side effects**

Undesirable side effects of this medical device, when the device is properly processed and used, are extremely rare. If you become aware of any adverse effects – even in doubtful cases – please inform us and describe the accompanying circumstances and symptoms as precisely as possible.

No risks in the event of an overdose are known of.

**Interactions**

Interaction of the product with other medical devices, medicinal products or other substances is not known.

**Storage**

Store C-Bond in a cool, dry place and protect from light. The Product should be stored between 0° C - 22°C.

**Minimum shelf life**

3 years

Do not use after the expiry date.

**Variants**

[REF]WP4300 1 x 5.0 ml liquid

[REF]WP4230 1 x 10.0 ml liquid

**Disposal**

Unused material can be cured and disposed of with household waste. Uncured material should be disposed of as methacrylate waste.

**All serious incidents related to the device shall be reported to the manufacturer and to the competent authority of the Member State in which the user and/or the patient is established.**

We know how

**dent a pharm**  
Produktionsgesellschaft mbH  
Schusterring 35  
25355 Barmstedt/Germany  
Fon: +(49) 4123/9225-0  
Fax: +(49) 4123/9225-48  
www.dent-a-pharm.de  
info@dent-a-pharm.de

**Willmann & Pein GmbH**  
Schusterring 35  
25355 Barmstedt/Germany  
Fon: +(49) 4123/9228-0  
Fax: +(49) 4123/9228-49  
www.wp-dental.de  
info@wp-dental.de

Made in Germany

26.03.2024  
2024-REV-012

H315 causes skin irritation  
H317 may cause allergic skin reactions  
H319 causes serious eye irritation  
H335 may cause respiratory irritation  
H411 toxic to aquatic organisms, with long-lasting effect

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