Securacem®

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Securacem is an X-ray opaque glass-ionomer cement that has been specifically developed as an underfilling material underneath composite and amalgam fillings. Securacem increases the biocompatibility of composites and improves the border seal. Adhesion to the dentin and enamel is by way of the reaction of the carboxylate groups of polyacrylic acid with the calcium or hydroxylapatite and collagen. This direct chemical bonding ensures that far more healthy dental material is saved because undercutting is not necessary.

Classification

DIN EN ISO 9917-1:2004-04 requirements are fulfilled

Preparing the cavities

A maximum of healthy dentine must be retained when preparing the cavities. The cavity must be shaped in such a manner that the cement layer is not thinner than 0.5 mm at any point.

Pulp protection

Rinse the prepared cavities thoroughly with water and dry with oil-free blown air. Use calcium hydroxide preparations as a pulp-capping for very deep cavities.

Preparing the dentin

Clean the freshly prepared dentin surface, rinse thoroughly with water and dry with blown air. Avoid dehydration.

Mixing Securacem

Shake powder well before use! Use one level dosing scoop for two full drop of liquid (without air bubbles; hold the bottle vertically downwards and apply light pressure to squeeze out one drop). First mix only half the amount of powder, then continue with the other half of the remaining powder until the required consistency is achieved. Use only a small area of the mixing block in order to minimize water evaporation.

Application of Securacem

The mixed material is lightly pressed into the cavity. Securacem is workable for 1'40" minutes (incl. mixing). The material should not be handled any more when hardening sets in. The hardening time is 3'00" - 3'40" minutes.

Instructions for use

Only cover the cavity base with Securacem if amalgam is to be used. The amalgam can be applied after 3'40" minutes. The dentin should be completely covered with Securacem when composites are used as filling material.

Treat the entire cavity surface (incl. the cement surface) with a dentin/enamel adhesive (e.g. C-Bond) according to the manufacturer's instructions. The composite filling is applied according to the manufacturer's instructions. We recommend the use of Competence HC for this purpose.

Side-effects

Side-effects with this medical product are extremely rare if it is properly prepared and applied. Should you encounter undesirable side-effects - even in case of doubt - then please contact us and give a precise description of the circumstances and symptoms.

Processing instructions

Ensure that the processing instructions and precautionary measures described in the supplied product information are always strictly observed. Should the patient suffer undesirable reactions or any other unwanted incidents arise even though these instructions were observed, then please inform us and give a precise description of the circumstances.

Storage

Shelf life: See label! Opened packs can be kept for 6 months provided they are kept dry and stored at room temperature.

WARNING! Do not allow non-bonded material to come into contact with the eyes. Should this nevertheless inadvertently happen, then rinse immediately with plenty of water.

Please note

Instructions must be immediately cleaned under flowing water after they have been used.

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