

LiWa Akti Bond

Description:

LiWa Akti Bond is an adhesive that allows to combine the modellation with LiWa resins, even when the first part of the modellation has already been light-cured. To do that, the polymerized modellation should be activated at the joint with a layer of LiWa Akti Bond.

LiWa Akti Bond is part of a wide range of materials: resins in different consistencies, preformed resin patterns, useful complementary materials for different areas of use as well as functional sets.

Properties:

- Liquid
- Easy to use
- Activates polymerized LiWa modellation to make them suitable for adding further LiWa resin
- Light-curing

Advantages:

- Time-saving
- Allows for a step by step modellation of LiWa, with intermediate curing of the finished parts

Indication:

Adhesive for a bond between polymerized LiWa parts and LiWa resin.

Storage conditions: dry and protected from light at 0 - 22 °C

Shelf life: 3 years

Presentation:

WP5082: LiWa Akti Bond, liquid, 10 ml bottle

